Celsite®

B | BRAUN

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MAGYAR P 97-100

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SLOVENSKÝ P 105-108

en: Description of the Device: Celsite® Access Ports are totally implantable access devices. The system consists of an access port, with a silicone septum, which is connected to a catheter. The catheter may be made of Silicone (Si), Polyurethane (PUR) or Polyamide (PA) depending on the reference. Access to the system is via percutaneous injection using a special non-coring needle. Both the port and the catheter are radio-opaque.

MR Information

General Information

According to IEC Standard 60601-2-33(2008), the scanner must be operated in Normal Operating Mode (defined as the mode of operation of the MR system in which none of the outputs have a value that causes physiological stress to the patient):

The whole body averaged specific absorption rate (SAR) must be ≤ 2.0 W/kg

The head SAR must be < 3.2 W/kg

Device Information

MR Conditional

The Celsite® Access Ports was determined to be MR-Conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-08. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.

Non-clinical testing demonstrated that the Celsite® Access Ports is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

Static magnetic field of 3-Tesla and 1.5-Tesla

Maximum spatial gradient magnetic field of 710 Gauss/cm or less

Maximum whole body averaged specific absorption rate (SAR) of 2.9 W/kg for 15 minutes of scanning

MRI-Related Heating

In non-clinical testing, the Celsite® Access Ports produced the maximum temperature rise during MRI performed for 15-min (i.e., per pulse sequence) in 3-Tesla (Excite, Software G3.0-052B, General Electric Healthcare, Milwaukee, WI) MR systems, as:

MRI Condition	MR System Reported, Maximum Whole Body Average SAR (W/kg)	Calorimetry Value (W/kg)	Highest Temperature Change	Time for MRI (per pulse sequence)
3-T / 128-MHz	2.9	2.7	2.2 °C	15 min

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Celsite® Access Ports. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

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INSTRUCTIONS FOR USE
Celsite® ACCESS PORTS
PLEASE READ CAREFULLY
SUMMARY

I	- INDICATIONS	Page 6
II	- CONTRAINDICATIONS	Page 6
III	- POTENTIAL COMPLICATIONS (Immediate and delayed)	Page 6
IV	- WARNINGS	Page 6
V	- SPECIAL ATTENTION	Page 6
	V-1 PATIENT	
	V-2 CATHETER POSITION	
	V-2-1 Venous catheters	
	V-2-2 Percutaneously placed Arterial catheters (Celsite® Interventional with Anthron® catheter)	
	V-3 DRUG INCOMPATIBILITIES	
	V-4 NEEDLES	
	V-4-1 Puncture of the Celsite® Access Port	
	V-4-2 Puncture of the vein	
	V-5 OCCLUSION	
	V-6 INFECTION	
VI	- IMPLANTATION TECHNIQUES	Page 6
	VI-1 CENTRAL VENOUS IMPLANTATION TECHNIQUES	
	VI-1-1 Percutaneous Technique	
	VI-1-2 Surgical Cut-down Technique	
	VI-1-3 References requiring modifications to the above techniques	
	VI-1-3-1 Over-the-Wire Technique	
	VI-1-3-2 Brachial Implantation Technique	
	VI-1-3-3 ECG Technique	
	VI-1-3-4 Braunule Technique	
	VI-1-3-5 Valved Catheter	
	VI-1-3-6 Pre-connected catheter and port	
	VI-2 ARTERIAL IMPLANTATION TECHNIQUES (FOR H.A.I.T.)	
	VI-2-1 Surgical Implantation Technique	
	VI-2-2 Percutaneous Implantation Technique	
	VI-2-2-1 Catheter tip in gastro-duodenal artery	
	VI-2-2-2 Catheter tip in hepatic artery	
	VI-3 PERITONEAL AND PLEURAL IMPLANTATION TECHNIQUES	
	VI-3-1 Peritoneal – Surgical Implantation Technique	
	VI-3-2 Peritoneal – Percutaneous Implantation Technique	
	VI-3-3 Pleural – Percutaneous Implantation Technique	
	VI-4 SPINAL IMPLANTATION TECHNIQUES	
VII	- USE OF THE PORT AND CATHETER	Page 8
	VII-1 VENOUS PORTS	
	VII-2 ARTERIAL PORTS	
VIII	- MAINTENANCE OF THE PORT AND CATHETER	Page 8
	VIII-1 RINSING AND HEPARINISATION	
	VIII-1-1 Venous ports	
	VIII-1-2 Arterial ports	
	VIII-1-3 Epidural or Intra-theal ports	
	VIII-1-4 Peritoneal ports	
	VIII-1-5 Pleural ports	
IX	- DURATION OF IMPLANTATION AND REMOVAL OF THE PORT	Page 8
	IX-1 DURATION OF IMPLANTATION OF THE SYSTEM	
	IX-2 REMOVAL OF THE SYSTEM	
	IX-3 SPECIAL ATTENTION: BABYPORT®, BABYPORT® S, BABYPORT® PC AND SMALL PORTS	
X	- HIGH PRESSURE INJECTION	Page 8

I – INDICATIONS

- Venous Access Ports are intended to be used for repeated intra-venous administration of, for example, chemotherapy, antibiotic and anti-viral drugs, parenteral nutrition, blood sampling or transfusion. Certain Celsite® Access Ports may be used for Contrast Enhanced Computerised Tomography (CECT) using high pressure injection (see 'X – High pressure injection'). For full list of validated references please refer to the tables on pages 109 and 110.
- Arterial Access Ports are intended to be used for intra-arterial administration of chemotherapy.
- Epidural or intra-theat ports are intended to be used for spinal administration of pain relieving drugs.
- Peritoneal ports are intended to be used for intra-peritoneal administration of chemotherapy.
- Peritoneal/Plural ports are intended to be used for intra-peritoneal administration of chemotherapy, hydration, drainage of malignant ascites, or drainage of malignant pleural effusion.

II – CONTRAINDICATIONS

The following are contraindications to port placement:

- Known infection, bacteraemia or septicaemia.
- Known allergy to any of the materials contained in the access port or catheter.
- If the medications to be used in the access port are incompatible with any of the materials contained in the access port or catheter.
- If the patient's anatomy does not allow the insertion of the catheter into the chosen access site or if the patient has had previous radiotherapy in the chosen area.
- Previous venous thrombosis.
- Heparin induced thrombocytopenia (Celsite® Interventional with Anthon® catheter only).

III – POTENTIAL COMPLICATIONS (Immediate and delayed)

Air embolism Angiospasm (arterial) Angitis (arterial) Aneurysm/false aneurysm (arterial) Brachial plexus injury Cardiac arrhythmia, tamponade, perforation Catheter disconnection, rupture or fragmentation Catheter occlusion	Cerebral infarction or thrombosis (arterial) Embolism Explantation of the port secondary to inflammatory reaction Fibrin sheath formation Haematoma Haemothorax Heparin Induced Thrombocytopenia (HIT) Infection/Sepsis	Liver failure (arterial) Pericarditis (arterial) Pneumothorax Port/Catheter migration Thromboembolism/Thrombophlebitis Thrombosis Ulcer related to cytotoxic drugs (arterial)
---	--	---

Before implantation of the access port system these, and other, well documented potential complications should be considered.

IV – WARNINGS

During storage avoid freezing or excessively high temperatures, storage at room temperature is recommended. Anthon® catheters should be protected from light. Celsite® Access Ports are for single use only, do not re-sterilise the product, or any of the components, and destroy after use. The device and its accessories are not re-usable nor designed to be re-used. Any re-use would definitely compromise the performance and safety of the device.

The product is sterilised using Ethylene Oxide and is sterile and non-pyrogenic in unopened, undamaged, individual packages.

Refer to unit package label for actual content list.

Check the expiry date before use and never implant a device past the expiry date.

Do not remove from the unit carton until ready for use and ensure that all handling and use is under strict aseptic conditions.

This product should not be implanted or accessed except by personnel experienced in the technical and clinical aspects of access ports.

Special Warning for Anthon® catheter:

- The Anthon® catheter is made of PUR with a heparin coating.
- The Anthon® catheter should not be used for blood sampling for coagulation studies as the results may be affected by the release, however small, of heparin.
- As for all heparin, regular supervision of platelet levels is required during the use of the catheter due to the risk of heparin induced thrombocytopenia (HIT).

V – SPECIAL ATTENTION

V-1 PATIENT

Ensure that use of the access port is indicated taking into consideration the patients' anatomy, needs and proposed use of the port.

Check the condition of the skin over the proposed implantation and puncture sites.

For venous catheter placement, place the patient slightly inclined head down.

V-2 CATHETER POSITION

To ensure that the system works correctly, there should be no kinking of the catheter.

It is recommended that the catheter is tunnelled sub-cutaneously to the port.

The position of the tip of the catheter may change when the patient moves.

V-2-1 Venous catheters

a) The catheter tip should always be positioned in the Superior Vena Cava at the entrance to the right atrium taking into account possible catheter tip movement (see Fig. 1, page 115). A radiograph should be taken to verify catheter position and exclude haemothorax or pneumothorax.

b) Particular attention should be paid when the catheter is to be implanted via the sub-clavian route. It is recommended that the catheter be inserted outside the costo-clavicular space (see Fig. 2, page 115). Catheter ruptures have been observed with the sub-clavian route, the associated risk of extravasation of the infused drugs and embolisation of the distal extremity may have serious consequences. This is due to rupture of the catheter secondary to the catheter being pinched in the costo-clavicular space (Pinch-Off syndrome). It is necessary to pay particular attention to catheters used for long periods or for ambulatory patients rather than short term catheters. Catheter rupture can occur with silicone and polyurethane catheters.

The following clinical signs may suggest catheter pinching:

- The patient needs to lift the arm to permit infusion.
- Intermittent malfunction of the catheter, such as difficulty with aspiration or infusion.
- Sub-clavicular pain or swelling during infusion.
- Palpitations or chest discomfort may indicate catheter fracture.

Remove any sub-clavian catheter which presents any of the signs of Pinch-Off.

A radiograph with contrast injection may be useful to detect catheter abnormalities at the costo-clavicular level.

V-2-2 Percutaneously placed Arterial catheters (Celsite® InterVenional with Anthon® catheter)

a) It is recommended that the catheter be placed via the femoral artery (see Fig. A, page 114). The entry point into the artery should be below the inguinal ligament. The access port should be positioned on the patient's thigh, medial to the iliac crest, to avoid kinking of the catheter (see Fig. B, page 114). The catheter should not traverse the inguinal ligament sub-cutaneously.

b) If needed, access may be via the axillary (see Fig. B, 2, page 114), acromioclavicular (see Fig. B, 3, page 114) or brachial artery (see Fig. B, 4, page 114). Care should be taken with these routes in patients with diabetes, hyperlipaemia and hyper-tension. The risk of cerebral embolism is higher in these patients.

c) Placement in the sub-clavian artery is not recommended as this may lead to catheter displacement and subsequent fracture, and carries an increased risk of cerebral thrombosis or infarct in case of catheter thrombosis.

V-3 DRUG INCOMPATIBILITIES

Never use incompatible drugs together or in sequence: the port should be rinsed of all residues with 10 mL sodium chloride (NaCl) 0.9% between each medication.

V-4 NEEDLES

V-4-1 Puncture of the Celsite® Access Port

Always use special bevelled needles, such as Surecan®, Winged Surecan® or Cytocan, which are adapted for use with implantable access port systems.

Normal hypodermic needles will damage the septum and may cause leakage of the system or blockage due to small silicone particles.

V-4-2 Puncture of the Vein

Certain references contain a Safety Seldinger needle. The Safety Seldinger needle is designed to reduce the risk of accidental needle stick injury. Care must be taken to avoid needle sticks in accordance with local Infection Control standards (and avoid the risk of exposure to contaminated blood).

Access the vein with the needle in the normal manner. After withdrawal of the Safety Seldinger needle activate the safety clip. To do this grasp the clear clip housing, which houses the metal safety clip, with your free hand. Advance the clear clip housing forward along the Safety Seldinger needle towards the needle level in one continuous motion. Advance the clip housing over the needle level activating the safety mechanism. The clip will remain on the needle level. The clip housing will separate from the needle and dispose. Dispose of the needle and the clear clip housing immediately into sharps container.

Do not attempt to re-insert the needle into the patient or to re-cap the needle with the protective guard.

V-5 OCCLUSION

Always use a syringe equal to or greater than 10 mL volume when rinsing the port to avoid generating excessive pressures in case of occlusion of the system.

In case of obstruction of the system never try to clear the blockage using a fluid under high pressure which carries the risk of catheter fracture and migration. According to validated protocols, and under medical supervision, 2 mL of 70% alcohol may be used to aid the unblocking of silicone catheters when the blockage is due to lipid deposits. The use of alcohol with PUR catheters is not recommended.

According to validated protocols, and under medical supervision, 2 mL of Hydrochloric acid (HCl) 0.1 mol/L may be used to aid the unblocking of both silicone and PUR catheters when the blockage is due to mineral deposits.

V-6 INFECTION

In case of infection, if appropriate treat with antibiotic drugs. If this fails, or is not appropriate, the catheter and port should be removed.

VI – IMPLANTATION TECHNIQUES

General considerations :

- To prevent shearing of the guide-wire never withdraw the J guide-wire through the Seldinger needle.
- Remove the guide-wire and the dilator together. Do not remove the guide-wire through the dilator as this may result in the guide-wire unravelling.
- To prevent shearing of the catheter, never withdraw the epidural/spinal catheters through the Tuohy needle.
- The recommended NaCl 0.9% (with or without heparin according to local protocols) is used to flush, and lock the port and catheter during implantation (see Fig. 3, page 115).
- During implantation ensure that the catheter is not damaged by unguarded forceps, suture needle or other sharp instruments.
- If using the Celsite® Concept ports (with silicone inserts) it is recommended that the suture be placed within the square indicated on the silicone.
- The catheter should be mounted on the exit cannula along its axis (see Fig. 11 a-b-c, page 115) and not across the axis and should be completely mounted on the exit cannula before the connection ring is slid over the catheter.

VI-1 CENTRAL VENOUS IMPLANTATION TECHNIQUES

VI-1-1 Percutaneous Technique

- Insert the Seldinger needle into the chosen vein; verify the position by observing blood reflux (see Fig. 4, page 115).
- Insert the J guide-wire into the vein, when the correct position is reached, withdraw the needle (see Fig. 5, page 115).
- Thread the assembled dilator and peelable sheath over the guide-wire using a twisting motion to pass through the skin planes (see Fig. 6, page 115).
- Remove the dilator and guide-wire from the vein and insert the catheter through the peelable sheath (see Fig. 7, page 115) to the desired position at the junction of the Superior Vena Cava and the right atrium (see Fig. 1, page 115). Check catheter tip position using fluoroscopy.
- Remove the peelable sheath from the vessel by pushing down on the 2 handles of the sheath, then peel the sheath from the catheter ensuring that the catheter is not dislodged from the vessel (see Fig. 8, page 115).
- Prepare the port pocket at the chosen site, the port should lie approximately 1/2 – 1 cm below the skin surface away from the injection site. Cut the luer connector from the catheter and attach to the tunnelling rod (see Fig. 9, page 115).
- Tunnel the catheter from the puncture point to the port pocket (see Fig. 10, page 115); ensure there is no kinking of the catheter. Cut the excess catheter (at right-angles) prior to connecting to the port.
- Catheter port/connection:

– Celsite®, Celsite® Concept, Celsite® Discret:

Slide the connection ring over the catheter, firmly push the catheter onto the exit cannula ensuring the catheter covers the length of the exit cannula, slide the connection ring over the catheter and exit cannula. The connection ring should be in contact with the port (see Fig. 11c, page 115).

– Celsite® IMPLANTOFIX:

Slide the screw connector over the catheter, firmly push the catheter onto the exit cannula ensuring the catheter covers the length of the exit cannula, slide the screw connector over the catheter and exit cannula. The screw connector should be tightened by hand and the spanner used to give the final tightening to the securing device, until an audible click is heard, to ensure correct attachment between the port and catheter.

– Celsite® Interventional with Anthon® catheter:

If necessary, cut the anti-kink portion of the connection ring to the desired length (a minimum of 1.5 cm should remain), then slide the connection ring over the catheter, firmly push the catheter onto the exit cannula ensuring the catheter covers the length of the exit cannula, slide the connection ring over the catheter and exit cannula. The connection ring should touch the port.

j) Confirm catheter connection by gently pulling on the catheter.

k) Insert the port into the port pocket, paying attention not to kink the catheter.

l) Confirm catheter patency (by ensuring that both aspiration and injection are possible), ensure that the skin incision is 1 cm from the injection site.

m) It is possible to fix the access port to the fascia with sutures.

– Celsite® DRAINAPORT:

Both the connection ring and the catheter must be threaded together on the exit cannula to facilitate connection.

VI-1-2 Surgical Cut-down Technique

- The skin is incised and the tissue dissected (see Fig. 12, page 115).
- The selected vessel is identified and tapes passed under it. The vessel is incised with a scalpel or with scissors (see Fig. 13, page 115).
- The previously flushed catheter is inserted into the vessel, using either a small vessel dilator or vein lifter (see Fig. 14, page 115).
- Follow from point e) – l) of the preceding paragraph: 'Percutaneous Technique'.

VI-1-3 References requiring modifications to the above techniques

VI-1-3-1 Over-the Wire Technique

This technique can only be used with the references ST301OTW and Celsite® IMPLANTOFIX references 04438604, 04438647.

Follow steps a) and b) above: 'Percutaneous Technique'.

c) Thread the dilator over the guide-wire using a twisting motion to pass through the skin planes. Remove the dilator (see Fig. 15, page 115).

Flush the catheter with sodium chloride (NaCl) 0.9 %, and pass the catheter over the guide wire to the desired position at the junction of the Superior Vena Cava and the right atrium (see Fig. 1, page 115).

e) Follow from point e) – l) above: 'Percutaneous Technique'.

VI-1-3-2 Brachial Implantation Technique

Brachial ports and catheters can be placed using either the OTW (Over-the Wire), Percutaneous or Surgical Cut-down Technique, (see Fig. E, page 114). Select the appropriate vein on the anterior aspect of the elbow, normally the Basilic or Cephalic vein. Ensure that the access port is not placed over the mid-line of the selected vein. The port and catheter should be positioned above the elbow on the internal surface of the arm.

VI-1-3-3 ECG Technique

This technique can only be used with the ECG references which have specially adapted accessories.

Percutaneous Technique: (see Fig. E, page 114)

- Place ECG pads on the patient's chest.
- Follow steps a) to c) above: "Percutaneous Technique".
- Note: A separate J-guide wire (0.035") is needed for access to the vein.
- Prepare the catheter by inserting the specially marked J-guide wire, which is contained in the kit, into the catheter.
- Remove the J-guide wire from the sheath and insert the catheter and guide-wire into the peelable sheath. The guide-wire should be pulled back into the catheter during placement into the vein.
- When approximately 10-13 cm of the catheter has been placed into the vein the guide-wire should be advanced into the catheter until the mark on the guide wire reaches the level of the luer connector on the catheter.
- When the mark on the guide-wire is at the level of the luer connector the tip of the guide-wire is in the desired position at the tip of the catheter.
- The sterile ECG cable is attached to the guide-wire, connected to the monitor via a universal adapter and switched to internal monitoring mode. The catheter and guide-wire are advanced using the ECG image to monitor progress.
- As the catheter approaches the right atrium the P wave starts to elevate. Maximal P wave height is reached and maintained when the catheter enters into the right atrium. Note that the amplitude of the P wave will decrease if the catheter is advanced down into the inferior vena cava. After identifying the area where the P wave begins to develop its maximal amplitude (which corresponds anatomically to the junction between superior vena cava and the right atrium) advance the catheter a further 2 cm. This is the final position of the catheter tip. As the procedure is performed with the patient in the supine position this catheter position allows for the 2-3 cm cranial movement of the catheter which occurs when the patient is upright.
- Remove the guide-wire while maintaining the catheter in place, flush and clamp the catheter.

Follow from point e) - l) above: "Percutaneous Technique".

A post-operative radiograph should be taken to confirm catheter position and exclude haemothorax or pneumothorax.

Surgical Cut-down Technique: (see Fig. F, page 114)

This technique combines the Surgical Cut-down Technique with the ECG technique described above.

If no ECG change is noted catheter tip position has to be confirmed with fluoroscopy as catheter tip malposition is to be suspected.

The Celsite® Access Port with ECG technique should not be placed if any of the following is applicable to the patient:

- The presence of a cardiac pacemaker.
- Known cardiac abnormalities, or heart rhythm disturbances.
- Access through the Inferior Vena Cava.

With regard to the ECG lead:

Safety regulations according to IEC 601-1 and NF EN 60601-2-25.

Protection against the risks of current leakage:

Only use Class 1 ECG monitor CE marked, with floating input (Type CF).

Connect additional potential equalisation (earth) to the monitor.

Protection against the risks of electrostatic discharge:

- Use only in rooms protected from electrostatic discharge (conductive flooring).
- Theatre shoes should be made of a recommended polymer material.
- Do not allow the J guide-wire to come into contact with external conductive parts.
- Before use, check that all connections are secure.
- Do not autoclave.
- Do not immerse in liquid.
- Repairs must be carried out only by the manufacturer.

VI-1-3-4 Braunule Technique

a) Puncture the vein with the Splittocan® needle percutaneously or through a skin incision.

b) Advance the catheter with the chlorhexidine (CHL) 0.9 %. Remove the stylette from the cannula, and insert the catheter through the cannula to the desired position at the junction of the Superior Vena Cava and the right atrium.

c) Gently withdraw the cannula over the catheter.

d) Follow from point e) - l) above: "Percutaneous Technique".

For the ST301M and ST305M, remove the mandrin from the catheter prior to cutting the Y connector off the catheter.

VI-1-3-5 Valved catheter

Do not place the catheter over the guide-wire, as this may cause damage to the valve.

The catheter must be tunneled from the puncture point to the port pocket; the tunneling rod must not be placed through the valve.

VI-1-3-6 Pre-connected catheter and port

The catheter should be measured and trimmed to the correct length before insertion in the vein.

The catheter must be tunneled from the port pocket to the puncture point.

VI-2 ARTERIAL IMPLANTATION TECHNIQUES (FOR H.A.I.T.)

H.A.I.T. = Hepatic Arterial Infusion Therapy

VI-2-1 Surgical Implantation Technique

- Perform a pre-operative arteriogram.
- Any accessory arteries to adjacent organs must be embolised.
- The catheter is sutured to the chosen artery using the suture beads at the tip of the catheter.
- Perfusion of the organ should be verified after placement of the catheter.
- Tunnel the catheter to the port pocket, the port should lie approximately 1/2 - 1 cm below the skin surface. Flush the port with sodium chloride (NaCl) 0.9 %.

VI-2-2 Percutaneous implantation Technique

Percutaneous placement of an arterial catheter is generally via the femoral artery. Access may also be via the axillary, acromioclavicular or brachial artery, and can be made either percutaneously or by surgical cut-down, depending on physician's choice.

An introducer should not be used during arterial procedures.

System implantation and catheter placement may be accomplished through a variety of techniques and the physician must choose the appropriate method for the patient and the intended result. Wider application, addition or alteration will be based on the experience of the physician. In the following chapter we present two of the most commonly used techniques for hepatic artery infusion.

VI-2-2-1 Catheter tip in gastro-duodenal artery (see Fig. C, page 114)

Note: The port/catheter reference R305-A55T has a side hole in the catheter 19 cm from the catheter tip. As the tip of the catheter will be embolised in place no flow occurs from the catheter tip. The taper of the catheter should be trimmed to the correct length (distal to the side hole) before placement in the gastro-duodenal artery so that the flow of medication is not obstructed and the side hole is in the proper hepatic artery.

A rigid 0.018" angiography guide-wire is required for placement of the tapered Anthon® catheter.

As the catheter is embolised in place, removal of the catheter percutaneously may not be possible.

a) An arteriogram should be performed to confirm the location of the tumour, the anatomy of the blood vessels and blood supply to adjacent organs.

b) Any accessory arteries to adjacent organs which arise distally to the proposed location of the side hole should be embolised.

c) Place an exchange guide-wire (minimum length 2 m, diameter 0.018") through the angiography catheter to the desired position in the gastro-duodenal artery.

d) Remove the angiography catheter and place the Anthon® catheter over the guide-wire (do not use an introducer).

e) Remove the guide-wire.

f) A check arteriogram should be performed to confirm catheter position.

g) The tip of the catheter should be fixed in place using embolisation coils, a gluing agent or other means of immobilising the catheter tip, and of embolising the gastro-duodenal artery. Access should be via a contra-lateral approach or through a micro-catheter which permits embolisation. Flow of medication is via the catheter side-hole.

h) In order to prevent displacement of the catheter tip during movement of the patient allow the catheter to sag or loop in the aorta.

Follow from point e) - l): "Percutaneous Technique".

The access port should be placed in the pelvic area ensuring that there is a solid base behind the port.

The port should not be placed over the abdominal cavity.

VI-2-2-2 Catheter tip in hepatic artery (see Fig. D, page 114)

a) An arteriogram should be performed to confirm the location of the tumour, the anatomy of the blood vessels and blood supply to adjacent organs.

b) Any accessory arteries to adjacent organs which arise distally to the proposed location of the catheter tip should be embolised.

c) Place a guide-wire through the angiography catheter to the desired position in the hepatic artery.

d) Remove the angiography catheter and place the Anthon® catheter over the guide-wire (do not use an introducer).

e) Remove the guide-wire.

f) A check arteriogram should be performed to confirm catheter position, function and adequate perfusion of the liver.

g) Ensure that a sufficient length of catheter is left in the aorta to prevent catheter tip displacement during patient movement.

Follow from point e) - l): "Percutaneous Technique".

The access port should be placed in the pelvic area ensuring that there is a solid base behind the port.

The port should not be placed over the abdominal cavity.

VI-3 PERITONEAL AND PLEURAL IMPLANTATION TECHNIQUES

The catheter and port may be implanted either surgically or percutaneously, and may be used for locoregional chemotherapy or for drainage of symptomatic malignant ascites.

VI-3-1 Peritoneal - Surgical implantation Technique

- During laparotomy the tip of the catheter is placed in the area of interest depending on the indication.
- The catheter may be sutured to the fascia if the cuffed catheter (Drainaport) is used, and then tunneled to the port pocket at the base of the ribs.
- Catheter function and position should be confirmed by flushing the port.
- The catheter is cut to the desired length and attached to the access port, taking care that the skin incision does not lie over the septum and that the catheter is not kinked anywhere along its length.
- Confirm catheter and port function before closing the skin.
- It is possible to fix the access port to the fascia with sutures.

VI-3-2 Peritoneal - Percutaneous implantation Technique

(Implantation accessories are available separately).

a) Identify the desired final position of the catheter tip using ultra-sound guidance.

b) With the 18G needle access the peritoneal cavity sub-cutaneously, taking care not to puncture the intestine. A curved needle may be used to aid entry into the peritoneal space.

c) Pass the 0.035" (0.89mm) guide-wire through the needle into the peritoneal cavity, remove the needle.

d) Pass the 12F dilator over the guide-wire to dilate the catheter route, followed by the 16F peelable sheath. Place the catheter into the desired position in the peritoneal space.

The use of a rigid catheter placed telescopically inside the peritoneal catheter and a super-stiff guidewire may facilitate its placement. Drain any fluid into a sterile bowl.

e) Create a pocket for the access port at the base of the ribs or over the iliac crest. If necessary further tunnel the catheter to the port pocket. It is advisable to have the cuff, if used, in the sub-cutaneous tissue so that tissue in-growth helps maintain catheter position.

f) Trim the catheter to the desired length and attach to the access port, take care that the skin incision does not lie over the septum and that the catheter is not kinked anywhere along its length.

g) Confirm catheter and port function before closing the skin.

h) It is possible to fix the access port to the fascia with sutures.

VI-3-3 Pleural - Percutaneous implantation technique

(Implantation accessories are available separately).

a) Identify the desired puncture point using ultra-sound (US) guidance. If the procedure is performed without US guidance, the puncture point will be marked by aspiration of pleural fluid through the needle used for injecting local anaesthetic.

b) With the 18G needle access the pleural space, taking care not to puncture the lung or liver.

c) Pass the 0.035" (0.89mm) guide-wire through the needle into the pleural cavity, remove the needle.

d) Pass the 12F dilator over the guide-wire to dilate the catheter route, followed by the 16F peelable sheath. Remove the guide-wire and dilator. Place the catheter into the desired position in the pleural space. Drain any fluid into a sterile bowl.

e) Create a subcutaneous pocket for the access port at the base of the ribs. Tunnel the catheter to the port pocket. Place the cuff in the sub-cutaneous tissue so that tissue in-growth helps maintain catheter position.

f) Trim the catheter to the desired length and attach to the access port. Take care that the skin incision does not lie over the septum and that the catheter is not kinked anywhere along its length.

g) Confirm catheter and port function before closing the skin.

h) It is possible to fix the access port to the fascia with sutures.

VI-4 SPINAL IMPLANTATION TECHNIQUES

a) Select the level for introduction of the catheter (L2-L3 or lower), and infiltrate the chosen area with local anaesthetic.

b) Select the position for the port pocket, this should be over a bony structure (the base of the ribs is usually chosen), and infiltrate local anaesthetic.

c) Make a small incision at the chosen entry site and enter the epidural or intra-thecal space using the Tuohy needle.

d) Introduce the catheter into the Tuohy needle using the insertion aid on the catheter, and advance into the chosen position, the catheter is marked at 5 cm intervals to aid placement. It is suggested that a minimum of 5 cm should be placed in the epidural space. The PU catheter should be advanced until the mandrin reaches the end of the needle, the mandrin should then be withdrawn while advancing the catheter into the intra-thecal or epidural space.

e) Remove the Tuohy needle from the puncture site while holding the catheter in a fixed position.

f) To verify the patency of the catheter attach the catheter to the catheter coupling device" (see Fig. 17, page 115). Attempt aspiration and inject 3-5 mL of dilute local anaesthetic for epidural placement. For spinal placement confirm the catheter position by ensuring aspiration of CSF (Cerebro-spinal-fluid) is possible.

g) Prepare the port pocket at the chosen site, the port should lie approximately 1/2 - 1 cm below the skin surface.

h) Thread the anti-kink device onto the catheter and attach the anti-kink device to the fascia with a suture.

i) Attach the catheter to the tunnelling rod, the route of the tunnelling rod may be anaesthetised with local anaesthetic prior to tunnelling. The catheter is then tunneled to the port pocket. It may be necessary to do this in two or more steps.

j) Cut the catheter, at right-angles, to the desired length allowing sufficient length for patient movement.

k) Slide the connection ring over the catheter; push the catheter onto the exit cannula of the chamber.

The connection ring should be tightened by hand and the spanner used to give the final tightening to the securement device.

l) The assembled port and catheter should again be checked for patency.

m) The chamber may then be sutured into the port pocket, taking care that the skin incision is not at the injection site.

* Attention: confirm that the catheter is pushed to the end of the catheter coupling device before closing to ensure correct function (see Fig. 17, page 115).

VII – USE OF THE PORT AND CATHETER

VII-1 VENOUS PORTS

Always verify that the port and catheter are functional by aspirating 2 mL of blood into a syringe and injecting 5 mL of sodium chloride (NaCl) 0.9% before attempting to start an infusion.

If aspiration of blood is not possible, attempt to slowly inject 2 mL of sodium chloride (NaCl) 0.9% into the port. If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected. In this case, device integrity should be verified using X-ray and contrast media studies.

During treatment, the infusion should be stopped immediately if any pain or swelling is noted or if blood return is absent.

If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected.

VII-2 ARTERIAL PORTS

Before attempting to start an infusion, verify that the port and catheter are functional. Using a syringe containing sodium chloride (NaCl) 0.9 %, attempt to inject into the port. Current practice is not to verify function by aspiration of blood into the syringe.

If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected.

VIII – MAINTENANCE OF THE PORT AND CATHETER

General recommendations:

- If heparinised sodium chloride (NaCl) 0.9% is used the system should be rinsed first with sodium chloride (NaCl) 0.9% alone, as some drugs may react with heparin and result in blockage of the port/catheter due to the formation of precipitates.

- Always rinse the system between administration of different drugs. Special care should be taken with drugs which have a higher risk of precipitation, with anti-coagulation agents, after blood sampling or after transfusion to reduce the risk of catheter occlusion.

- Failure to maintain the system may result in occlusion of the catheter.

VIII-1 RINSING AND HEPARINISATION

VIII-1-1 Venous ports

Rinsing of the access port is essential.

Rinse the port with 10 mL of sodium chloride (NaCl) 0.9% after each use and every 4 weeks when no treatment is being given.

If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected.

VIII-1-2 Arterial ports

Arterial ports and catheters should be rinsed with 10 mL of sodium chloride (NaCl) 0.9% followed by heparinised sodium chloride (NaCl) 0.9% after each treatment and regularly, according to clinical needs (eg every 4 weeks) when no treatment is being given.

If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected.

VIII-1-3 Epidural or Intra-thecal ports

Sodium chloride (NaCl) 0.9%, 0,5 mL may be used to flush the system after use. Heparin should not be used in an epidural or spinal port.

If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected.

VIII-1-4 Peritoneal ports

Peritoneal ports may be rinsed with 20 mL sodium chloride (NaCl) 0.9% after each treatment.

If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected.

VIII-1-5 Pleural ports

Pleural ports should be rinsed with 20 mL of heparinised sodium chloride (NaCl) 0.9%, initially this should occur weekly. The interval between rinsing may be extended according to clinical needs.

If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected.

IX – DURATION OF IMPLANTATION AND REMOVAL OF THE PORT

IX-1 DURATION OF IMPLANTATION OF THE SYSTEM

The port and catheter system should be removed at the end of treatment.

The duration of the system is mostly dependant on the catheter. Fracture is the most common complication of ageing. Risk of fracture may increase over time depending on the material of the catheter, the pathway, and the entry point of the catheter into the vessel.

IX-2 REMOVAL OF THE SYSTEM

When removing the system, care must be taken not to fracture the catheter. If the catheter is sutured into the vessel, the sutures should first be removed. Control the catheter while removing the system from the port pocket. Be vigilant if there is excessive resistance to the removal of the catheter. Catheters may become encapsulated and attached to the vein wall. Should this occur and it is not possible to remove the catheter without risking catheter fracture, or if the catheter is fractured, the advice of an Interventional Radiologist, Surgeon or other Physician with endo-luminal experience should be sought.

Arterial catheters - If the catheter has been immobilised in an artery, using either coils or gluing techniques, removal of the catheter may not be possible. Attempts at removal may result in fracture of the catheter within the artery and/or vessel injury.

IX-3 SPECIAL ATTENTION: BABYPORT®, BABYPORT® S, BABYPORT® PC AND SMALL PORTS

Following implantation of a port in a child particular attention should be given to the position of the catheter tip. Over time, as the child grows, the position of the catheter tip will change and move higher in the Superior Vena Cava. As for any catheter, a high tip position in the SVC predisposes the catheter to fibrin sheath formation and other mechanical difficulties.

The position of the catheter tip should be checked radiographically a minimum of every 12 months or more frequently if the child is growing quickly. The system should be changed when the tip reaches the level of T4.

X – HIGH PRESSURE INJECTION

SPECIAL ATTENTION

- Always verify that the CelSite® reference implanted is included in the table on page 109 and 110, only these references may be used for high flow rate/high pressure injection.

- Always verify that the port and catheter are functional by aspirating 2 mL of blood into a syringe and injecting 5 mL of sodium chloride (NaCl) 0.9% into the port/catheter before attempting to start an infusion of medication.

- Do not exceed the recommended pressure (325 psi- 22.4 bars) and flow rate as Access Port system failure may occur.

- Contrast media should be warmed to 37°C (98.6°F) before use. According to drug manufacturers recommendations. Failure to follow this recommendation will result in up to 50% lower flow rates and/or Access Port or injection system failure.

- Use only Winged Surecan® (without Y-site) or Angled Surecan® needles.

- Do not use needles which may not withstand high pressure.

- Use only 20G or 22G needles for Babyport® and Brachial ports. Use of 19G needles may lead to leakage of contrast media.

- Ensure that the needle is correctly placed in the port, securely taped to the skin and covered with an adhesive dressing before commencing high pressure injection.

- Catheter lengths more than 20 cm will result in reduction in flow rates.

- Depending on the technical characteristics of the injector system, the target flow rate might not be attained.

- The Access Port system should be flushed with 10 mL of sodium chloride (NaCl) 0.9% before, and after, using the port for CECT, followed by usual rinsing procedures.

		en: Recommended maximum flow rates (mL/s) ⁽¹⁾						Recommended maximum pressure setting (CT function) ⁽⁸⁾
		Contrast media at 37°C (98,6°F) ⁽³⁾						
		Viscosity 5,8 mPa.s (cP)** ⁽⁴⁾			Viscosity 11,4 mPa.s (cP)** ⁽⁵⁾			
		Winged / Angled Surecan® needle ⁽⁷⁾			Winged / Angled Surecan® needle ⁽⁷⁾			
		22G	20G	19G	22G	20G	19G	
	Babyport® - Babyport® PC (1,5 x 0,8 mm)	2	4		1	3		325 psi
	Brachial (1,6 x 1,1 mm)	2	4		1	3		
	Brachial L - Brachial R (1,6 x 1,1 mm)	2	4		1	3		
	Babyport® S (2,0 x 1,2 mm)	2	4		2	4		
Double port ⁽⁹⁾	ST405L (3,2 x 1,2 mm)	2	5	6	2	4	6	
Small Size Port ⁽¹⁰⁾	STL205P - STR205P (2,1 x 1,4 mm)	2	4	6	2	3	5	
	ST305P (2,1 x 1,4 mm)	2	4	6	2	3	4	
	ST305C - CR305-A5 (1,7 x 1,1 mm)	2	4	5	1	3	4	
	CR305-A6 (2,0 x 1,2 mm)	2	4	7	2	3	4	
	T/ST305 - T/ST205 - ST505 ST315 - ST215 - ST205F ECG (2,2 x 1 mm)	2	4	5	2	3	4	
	STL205F - STR205 F (2,2 x 1 mm)	2	4	5	2	3	4	
	ST305L - ST505L - ST205ECG - ST315L (2,8 x 1,1 mm)	2	4	5	2	3	5	
	ST305M (2,0 x 1,2 mm)	2	4	6	2	3	4	
	ST305H - ST505H (2,8 x 1,6 mm)	2	5	7	2	4	6	
Double port ⁽⁹⁾	ST401L (3,2 x 1,2 mm)	2	5	7	2	4	6	
Standard Size Port ⁽¹¹⁾	ST301M (2,0 x 1,2 mm)	2	5	8	2	4	5	
	ST301C - ST501C - ST201C - ST301OTW CR301-A5 (1,7 x 1,1 mm)	2	5	6	2	4	5	
	CR301-A6 (2,0 x 1,2 mm)	2	5	6	2	4	6	
	T/ST301F - ST311F T/ST201F - T/ST501F - ST201F ECG (2,2 x 1 mm)	2	5	6	2	4	6	
	T/ST301P - ST201P (2,1 x 1,4 mm)	2	5	6	2	4	6	
	T/ST301 - ST311 T/ST201 - T/ST501 - ST201ECG (2,8 x 1,1 mm)	2	5	6	2	4	6	
	STL201L - STR201L (2,8 x 1,1 mm)	2	5	6	2	4	6	
	ST201H - T/ST301H - ST311H - ST501H (2,8 x 1,6 mm)	2	5	7	2	5	7	
	STL201H - STR201H (2,8 x 1,6 mm)	2	5	7	2	5	7	
	ST301G - ST201G - ST501G (3,2 x 1,6 mm)	2	5	8	2	5	7	

* With a catheter of 20 cm⁽²⁾

** Or any contrast media with a similar viscosity at 37°C⁽⁶⁾

en: FOR COUNTRY UNDER CE MARK ONLY / fr: POUR LES PAYS SOUS MARQUAGE CE UNIQUEMENT / de: NUR FÜR LÄNDER MIT GÜLTIGKEIT DES CE ZEICHENS / es: SÓLO APLICABLE A LOS PAÍSES BAJO LA MARCA CE / it: SOLO PER PAESI IN CUI È NECESSARIO IL MARCHIO CE / sv: ENDAST FÖR LÄNDER MED CE-MÄRKNING / pt: APENAS PARA PAÍSES SUJEITOS A MARCA CE / zh: 仅限CE标志适用国家 / ru: ТОЛЬКО ДЛЯ СТРАН ПОД МАРКОЙ CE / nl: VOOR LANDEN MET CE-MARKERING / no: KUN FOR LAND MED CE MERKING / da: KUN FOR CE MÆRKEDE LANDE / fi: MAAT JOILLA ON CE-MERKINTÄ / pl: WYŁĄCZNIE DLA KRAJÓW UŻYWAJĄCYCH OZNACZENIE CE / el: ΜΟΝΟ ΓΙΑ ΧΩΡΕΣ ΟΠΟΥ ΑΠΑΙΤΕΙΤΑΙ ΠΙΣΤΟΠΟΙΗΣΗ CE / tr: SADECE CE BELGESİNİN GEÇERLİ OLDUĞU ÜLKELER İÇİN / cs: POUZE PRO ZEMĚ S OZNAČENÍM CE / ko: 유럽연합 통합규격마크 인증 받은 나라만 해당 / vi: CHỈ ÁP DỤNG CHO NHỮNG NƯỚC ĐẠT CHỨNG NHẬN EC / ro: DOAR PENTRU ȚĂRILE CU MARCAJ CE / bg: ЗА ДЪРЖАВИ, ПРИЛАГАЩИ СЕ МАРКИРОВКА / et: CE MÄRK-GISTUS / hu: ORSZÁG SZERINTI CE ENGEDÉLY / lt: CE ŠALIMS / sk: IBA PRE KRAJINU S OZNAČENÍM CE

		Recommended maximum flow rates (mL/s) ⁽¹⁾						Recommended maximum pressure setting (CT function) ⁽⁸⁾
		Contrast media at 37°C (98,6°F) ⁽³⁾						
		Viscosity 5.8 mPa.s (cP)** ⁽⁴⁾			Viscosity 11.4 mPa.s (cP)** ⁽⁵⁾			
		Surecan® Safety II			Surecan® Safety II			325 psi
		22G	20G	19G	22G	20G	19G	
Double port ⁽⁹⁾	Babyport® - Babyport®PC (1,5 x 0,8mm)	2	4	-	1	3	-	325 psi
	Brachial (1,6 x 1,1mm)	2	5	-	1	4	-	
	Brachial L - Brachial R (1,6 x 1,1mm)	2	5	-	1	4	-	
	Babyport® S (2,0 x 1,2mm)	2	5	-	1	4	-	
	ST405L (3,2 x 1,2mm)	2	5	8	1	4	6	
	STL205P - STR205P (2,1 x 1,4mm)	2	5	8	1	4	5	
	ST305P (2,1 x 1,4mm)	2	5	8	1	4	5	
	ST305C - CR305-A5 (1,7 x 1,1mm)	2	5	7	1	3	5	
	CR305-A6 (2,0 x 1,2 mm)	2	5	7	1	3	5	
T/ST305 - T/ST205 - ST505 ST315 - ST215 - ST205F ECG (2,2 x 1mm)	2	5	8	1	4	6		
STL205F - STR205F (2,2 x 1mm)	2	5	8	1	4	6		
ST305L - ST505L - ST205ECG - ST315L (2,8 x 1,1mm)	2	5	8	1	3	6		
ST305M (2,0x 1,2mm)	2	5	8	1	3	6		
ST305H - ST505H (2,8 x 1,6mm)	2	6	9	1	4	6		
ST401L (3,2 x 1,2mm)	2	5	8	1	4	6		
Standard Size port ⁽¹¹⁾	ST301M (2,0 x 1,2mm)	2	6	7	1	4	6	
	ST301C - ST501C - ST201C - ST301OTW CR301-A5 (1,7 x 1,1mm)	2	5	6	1	3	5	
	CR301-A6 (2,0 x 1,2 mm)	2	5	6	1	3	5	
	T/ST301F - ST311F T/ST201F - T/ST501F - ST201F ECG (2,2 x 1mm)	2	6	7	1	4	5	
	T/ST301P - ST201P (2,1 x 1,4mm)	2	5	7	1	4	6	
	T/ST301 - ST311 T/ST201 - T/ST501 - ST201ECG (2,8 x 1,1mm)	2	6	7	1	4	6	
	STL201L - STR201L (2,8 x 1,1mm)	2	6	7	1	4	6	
	ST201H - T/ST301H - ST311H - ST501H (2,8 x 1,6mm)	2	6	8	1	4	7	
	STL201H - STR201H (2,8 x 1,6mm)	2	6	8	1	4	7	
ST301G - ST201G - ST501G (3,2 x 1,6mm)	2	6	9	1	4	6		

⁽¹⁾With a catheter of 20 cm⁽²⁾

⁽²⁾Or any contrast media with a similar viscosity at 37°C⁽⁸⁾

en: Translations : see page 111 / fr: Traductions : Voir page 111 / de: siehe Seite 111 / es: Traducción: ver página 111 / it: Traduzioni: vedere a pagina 111 / sv: Översättningar: Se sidan 111 / pt: Traduções: consulte a página 111 / zh: 翻译: 请参阅第 111 页 / ru: Перевод: см. на стр.111 / nl: Vertalingen: zie pagina 111 / no: Oversettelser: se side 111 / da: Oversættelse: se side 111 / fi: Käännös: katso sivu 111 / pl: Tłumaczenia: patrz strona 111 / el: Μεταφράσεις: Βλέπε σελίδα 111 / tr: Çeviriler : Sayfa 111'e bakınız / cs: Překlady: viz str. 111 / ko: 번역: 111쪽 / vi: Phần dịch: xem trang 111 / ro: Pentru traducere consultați pagina 111 / bg: Превод: моля виж стр. 111 / et: Tõlked : vt. Lk. 111 / hu: fordítás: lásd a 111. oldalon / lt: vertimų žiūrėkite psl. Nr.111 / sk: Preklady: viď str. 111.

NEW

Celsite® ST205L / ST205P / ST205H CT Marking of Epoxy Access Ports

Special Attention: FOR COUNTRY UNDER CE MARK ONLY

The following venous Access Ports are marked with "CT" to indicate compatibility with power injection. This marking is visible on X-ray.

All characteristics of these devices are summarized below :

Reference	Port/catheter material	Dead volume of port	Dead volume of catheter/cm	Implantation route	Possible implantation techniques
ST205L	Epoxy Resin / Silicone	0.25 ml	0.02 ml	Venous	Surgical, Seldinger
ST205P	Epoxy Resin / PUR	0.25 ml	0.015 ml	Venous	Surgical, Seldinger
ST205H	Epoxy Resin / PUR	0.25 ml	0.02 ml	Venous	Surgical, Seldinger

		Recommended maximum flow rates (mL/s)*											
		Contrast media at 37°C (98,6°F)											
		Viscosity 5.8 mPa.s (cP)**			Viscosity 11.4 mPa.s (cP)**			Viscosity 5.8 mPa.s (cP)**			Viscosity 11.4 mPa.s (cP)**		
		Surecan® Safety II			Surecan® Safety II			Winged/Angled Surecan			Winged/Angled Surecan		
		22G	20G	19G	22G	20G	19G	22G	20G	19G	22G	20G	19G
Small Size port (10)	ST205P (2,1 x 1,4mm)	2	5	8	1	4	5	2	4	8	2	3	4
	ST205L (2,8 x 1,1mm)	2	5	8	1	3	6	2	4	5	2	3	5
	ST205H (2,8 x 1,6mm)	2	6	9	1	4	6	2	5	7	2	4	6
*With a catheter of 20 cm (2)													
**Or any contrast media with a similar viscosity at 37°C (6)													

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en: All characteristics of these devices are summarized below :

Reference	Port/catheter material	Dead volume of port	Dead volume of catheter/cm	Implantation route	Possible implantation techniques
BT301P	Polysulphone / PUR	0.50 ml	0.015 ml	Venous	Surgical
BT305P	Polysulphone / PUR	0.25 ml	0.015 ml	Venous	Surgical

	Recommended maximum flow rates (mL/s)* (1)											
	Contrast media at 37°C (98,6°F) (3)											
	Viscosity 5.8 mPa.s (cP)** (4)			Viscosity 11.4 mPa.s (cP)** (5)			Viscosity 5.8 mPa.s (cP)** (4)			Viscosity 11.4 mPa.s (cP)** (5)		
	Surecan® Safety II			Surecan® Safety II			Winged/Angled Surecan (7)			Winged/Angled Surecan (7)		
	22G	20G	19G	22G	20G	19G	22G	20G	19G	22G	20G	19G
BT301P (2,1 x 1,4 mm)	2	5	7	1	4	6	2	5	6	2	4	6
BT305P (2,1 x 1,4 mm)	2	5	8	1	4	5	2	4	6	2	3	4

*With a catheter of 20 cm (2)
**Or any contrast media with a similar viscosity at 37°C (6)

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en: Special Attention

CT Marking of Epoxy Access Ports

Contrast Enhanced Computerized Tomography (CECT) – refer to I – Indications

The following venous access ports are marked with “CT” to indicate compatibility with power injection. This marking is visible on X-ray. Please refer to page 109 for recommended maximum flow rates (mL/sec).

Paediatric	Small / Double	Standard / Double
Babyport®	STL205P, STR205P	ST201C
Babyport® PC	T/ST205, STL205F, STR205F	T/ST201F
Brachial	ST215	ST201F ECG
Brachial L, Brachial R	ST205F ECG	T/ST201, STL201L, STR201L
Babyport® S	ST205ECG	ST201ECG
	ST405L	ST201H, STL201H, STR201H
		ST201G
		ST201P
		ST401L

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