AESCULAP®

en Instructions for use/Technical description

Neuro-Patch®

USA Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusaifus.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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1. About this document

Note

General risk factors associated with surgical procedures are not described in these instructions for use.

1.1 Scope

These instructions for use apply for the following products:

Art. no.	Designation	Content
1064002	Neuro-Patch 12 cm x 14 cm	1 pc
1064010	Neuro-Patch 6 cm x 14 cm	1 pc
1064020	Neuro-Patch 8 cm x 9 cm	1 pc
1064029	Neuro-Patch 6 cm x 8 cm	1 pc
1064037	Neuro-Patch 4 cm x 10 cm	1 pc
1064040	Neuro-Patch 5 cm x 6 cm	1 pc
1064045	Neuro-Patch 4 cm x 5 cm	2 pcs
1064053	Neuro-Patch 2 cm x 10 cm	2 pcs
1064061	Neuro-Patch 1.5 cm x 3 cm	2 pcs
1064110	Neuro-Patch 4 cm x 5 cm	1 pc
1064122	Neuro-Patch 2 cm x 10 cm	1 pc
1064123	Neuro-Patch 1.5 cm x 3 cm	1 pc

2. Clinical use

Note

The short summary of safety and clinical performance of the product (SSCP) is available in the European Database for Medical Products (EUDAMED).

As long as EUDAMED is not available the SSCP can be downloaded at <http://www.bbraun.com/neurosurgicalpassport>

2.1 Product description

Neuro-Patch is a non-absorbable, fine-fibrillar, microporous nonwoven fabric made from high-purity, aliphatic poly-ester urethane. The structure of Neuro-Patch is characterized by open micropores on the surface of the nonwoven fabric, which are aimed at significantly facilitating the immigration of fibroblasts.

2.2 Areas of use and limitations of use

2.2.1 Intended use

Neuro-Patch is used in neurosurgery as dura mater replacement.

2.2.2 Indications

Note

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

Usage:

- For covering cerebral and cerebellar dura defects
- For cerebral decompression surgery when there is elevated intracranial pressure
- For covering spinal dura defects
- For spinal decompression surgery

2.2.3 Contraindications

Do not use:

- In infected areas
- In open cerebrocranial traumata
- In open spina bifida
- In case of known hypersensitivity against implant materials; for fixation materials please note the corresponding instructions for use
- In any application area that is not mentioned in "Indications"

2.3 Risks, adverse effects and interactions

Potential complications that the manufacturer is currently aware of:

- Infection
- Cerebrospinal fluid leakage
- Adhesion
- Foreign body reaction

Note

The points mentioned above include potential clinical consequences.

No risks, side effects and interactions as a result of comorbidities of the patient have been identified.

2.4 Safety information

2.4.1 Clinical user

General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Always follow the safety advice and information given in the instructions for use.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Keep instructions for use accessible to the user.

Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

Aesculap is not responsible for complications caused by:

- incorrect indication or implant selection
- incorrect surgical technique
- incorrect combination of implant components
- combination with components of other manufacturers not approved by Aesculap
- exceeding the limitations of the treatment method or non-observance of essential medical precautions

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

2.4.2 Sterility and storage

The product is EO sterilized and wrapped in sterile packaging.

- ▶ Store implant components in their original packaging. Remove them from their original protective packaging only just prior to implantation.
- ▶ Do not use products from open or damaged sterile packaging.
- ▶ Do not use the product after its use-by date.
- ▶ Store the product at 25 ± 5°C.
- ▶ Do not reuse the product.

The reprocessing of the product affects its functionality. Risk of injury, illness or death due to soiling and/or impaired functionality of the product.

- ▶ Do not reprocess the product.

2.5 Patient education

Within the framework of the patient education, the relevant circumstances needed for consent must be explained to the patient in accordance with their level of understanding, pre-existing knowledge and need for information. This includes:

- Diagnosis, procedure and risk clarification
- Operative procedure
- Advantages and disadvantages of the procedure
- All alternative procedures that can be considered

2.6 Application

2.6.1 Documentation

The user shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of implants

Important information concerning the implanted product and the surgical procedure can be noted down on the implant card. The implant card is enclosed with the product.

Note

A patient information brochure and a guide on how to fill out the implant card are available online. The corresponding website is indicated on the implant card.

- ▶ Fill out the implant card and hand it to the patient together with the patient information brochure after the surgical procedure.

Each package contains additional labels showing the designation, article and lot number of the product.

- ▶ Use these labels for documentation in the patient's file.

The surgical procedure has been explained to the patient, and the patient's consent has been documented.

2.6.2 Implantation

- ▶ Choose the suitable implant size to fit the application area.
- ▶ Cut Neuro-Patch into the required shape.
 - ▶ The implant should be cut as closely as possible to the defect size to achieve tension-free embedding.
- ▶ Fix the Neuro-Patch using non-absorbable suture material (polyester, polypropylene).

Note

The use of atraumatic round body needles permits suturing without major damage to the implant.

In addition, fibrin glue can be used to achieve a sealing effect.

Note

During the application of Neuro-Patch in combination with bone cement, chemical damage of the patch material, depending on the application situation, cannot be excluded.

3. Disposal

Contaminated products pose a risk of infection!

- ▶ Observe national regulations when disposing or recycling the product, its components and their packaging.

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