



ASKINA[®] SILNET

SOFT SILICONE WOUND CONTACT LAYER

B. Braun Wound Care. Right. From the start.

Askina® SilNet

Soft silicone wound contact layer



Non - adherent: protects fragile granulation tissue

Atraumatic dressing change

Soft and conformable

Allows passage of exudate through the dressing

Askina® SilNet is made of a conformable non-woven material, coated on both sides with soft silicone.


It is used as a wound contact layer, it conforms to the wound surface and adheres safely to the surrounding skin.

The wound site is protected from mechanical disruption during dressing changes. The use of Askina® SilNet thus minimises the trauma associated with dressing change.

TIP FOR APPLICATION

Askina® SilNet can be used in combination with secondary dressings

- Askina® Foam
- Askina® Heel
- Askina® Pad
- As an interface under NPWT or to contain Prontosan® Gel/Gel X in a wound



INDICATIONS

- Skin tears
- Traumatic wounds
- Superficial partial thickness and deep partial thickness burns
- Fixation of skin grafts
- Donor sites
- Lacerations
- Skin Grafts

HOW TO USE

1. Remove the transparent liner.
2. Apply Askina® SilNet on the wound surface and gently remove the white liner.
3. Apply suitable absorbing secondary dressing, like Askina® Foam.
4. Its porous structure allows vertical passage of exudate into the secondary absorbing dressing, with no risk of maceration.

Soft silicone adhesive

Askina® SilNet is made of soft silicone, which provides safe and secure fixation, without stripping of epidermal cells during removal, unlike traditional acrylic adhesive.



B. Braun Ordering Information

Askina® SilNet

Size	Pcs/Pack	Article number
5 cm x 7.5 cm	10	5195710
10 cm x 7.5 cm	10	5197510
10 cm x 18 cm	10	5191810
20 cm x 30 cm	5	5192305

B. Braun Australia Pty Ltd | Bella Vista NSW | Tel. 1800 251 705 | info@bbraun.com | www.bbraun.com.au

B. Braun New Zealand | Customer Care 0800 227 286 | Fax (09) 373 5601 | www.bbraun.co.nz

BAUS OM E0562 03/20 | ©2020 B. Braun Australia