

# Drug Coated Balloon Technology within Coronary Angioplasty

BASIC TRAINING FOR  
CATH LAB STAFF

**B|BRAUN**  
SHARING EXPERTISE



OUR VISION

We  
**PROTECT**  
and  
**IMPROVE**  
the  
**HEALTH**  
of people  
around the world.

## CONTENTS

Evolution in Angioplasty

Challenges with Drug Eluting Stents

Case Studies

ISR

Small Vessel Disease

Bifurcations

DCB-only Approach

Late Lumen Enlargement

Consensus Guidelines

The Angioplasty Process

SeQuent® Please NEO

Handling Tips & Tricks

Advantages

What is the Key Message?

SeQuent® Please NEO Specification


Ordering Information

Technical Data

References

# EVOLUTION IN ANGIOPLASTY

It is important to see where technologies have come from in order to truly understand them.



# EVOLUTION IN ANGIOPLASTY

PCI is one of the most common medical procedures worldwide. PCI celebrated its 40th year anniversary in 2017 and since 1977, there has been a significant improvement in treatment technologies.

## Plain Old Balloon Angioplasty (POBA)

Although revolutionary in 1977, POBA produced unpredictable procedural outcomes due to vessel dissection, vessel recoil and high restenosis rates. The patient would often require repeat dilation or bypass surgery.

## Bare Metal Stents (BMS)

BMS produced more stable results with lower restenosis rates. However, challenges with this technology included early stent thrombosis and neointimal hyperplasia – proliferation of smooth muscle cells within the intima layer of the artery wall leading to a decrease in lumen diameter and vessel re-narrowing.

## Drug Eluting Stents (DES)


The introduction of DES significantly improved patient outcomes by delivering an anti-proliferative agent to the vessel wall. This lowered the rate of in-stent restenosis (ISR) and allowed the use of PCI in more advanced and complex disease. Whilst first generation DES were correlated to an increase in late stent thrombosis and patient death, second generation DES have been proven to be safe and effective<sup>3</sup>.

# CHALLENGES WITH DRUG ELUTING STENTS (DES)


## IN-STENT RESTENOSIS (ISR)

The re-narrowing of a blood vessel which was previously treated with a stent during PCI leading to restricted blood flow.

ISR remains a challenge in the interventional cardiology community, on average ISR occurs in 5% of all PCI cases<sup>5</sup>.



At time of PCI with BMS




6 weeks post intervention



The European Society of Cardiology have given DCBs a Class 1A recommendation for ISR. Such recommendations are based around the vast amount of supporting clinical data in randomised controlled trials<sup>7</sup>.

## ISR DCB-ONLY CASE STUDY



Pre-treatment


Post-treatment

4-month follow-up

- Patient:** Male, 55 years
- Indication:** ISR of BMS (3.5 x 15 mm) implanted 2 years previously
- Procedure:** **LESION PREPARATION**  
Pre-dilation PTCA balloon (3.5 x 15 mm)
- LESION THERAPY**  
DCB-only SeQuent® Please  
(3.5 x 20 mm) proximal lesion  
DCB-only SeQuent® Please  
(3.5 x 15 mm) distal lesion
- Post-treatment:** No dissection proximal or distal to the BMS, no residual stenosis present
- Follow-up:** A very successful result without the need for additional metal within the artery

Full case studies available on request.

## SMALL VESSEL DISEASE DCB-ONLY CASE STUDY



Pre-treatment


Post-treatment

5-month follow-up

- Patient:** Male, 74 years
- Pre-treatment:** Long stenosis of distal right coronary artery
- Procedure:** **LESION PREPARATION**  
Pre-dilation PTCA balloon (2.5 x 20 mm)  
Two inflations, distal and proximal of posterior descending
- LESION THERAPY**  
DCB-Only SeQuent® Please (2.5 x 20 mm x 2)  
distal and posterior descending
- Post-treatment:** No relevant stenosis or dissections present, good TIMI flow
- Follow-up:** No stenosis after 5 months and a very acceptable result using the DCB-only approach

Full case studies available on request.

# BIFURCATIONS DCB-ONLY CASE STUDY



Pre-treatment

Post-treatment

4-month follow-up

**Patient:** Male, 54 years

**Pre-treatment:** Stenosis of mid circumflex artery (CX) and its posterolateral branch (PL-CX)

**Procedure:** **LESION PREPARATION**

Pre-dilation PTCA Balloon (2.5 x 20 mm)

**LESION THERAPY**

DCB-Only SeQuent® Please  
(3.0 x 15 mm) of PL-CX  
(3.0 x 20 mm) of CX

**Post-treatment:** No residual stenosis present and good TIMI flow

**Follow-up:** A successful result without the need for any metal work within the artery

Full case studies available on request.

# DCB-ONLY APPROACH

## GO IMPLANT FREE

The DCB-only approach is the use of a DCB in place of a stent for a primary coronary intervention and therefore is the concept of 'going implant free' i.e. leaving no metal work behind within the artery.

## ADVANTAGES



Avoidance of stent related complications  
e.g. ISR or stent thrombosis




DAPT can be reduced to one month – important  
for HBR patients or those requiring surgery  
shortly after intervention<sup>8,9,10</sup>




Enables late lumen enlargement of the target  
vessel and therefore positive remodelling<sup>8</sup>

## WHAT IS LATE LUMEN ENLARGEMENT?




Following treatment using the DCB-only approach, the target vessel lumen diameter can in fact increase in size over time. This is because there is nothing implanted or left behind in the artery, allowing natural vasomotion to continue. Further, the vessel is not restricted in terms of diameter as it is with a stent.

Once a stent is implanted into the vessel, the vessel diameter can never increase, it is restricted by the metal stent struts and therefore can only ever decrease in size. Using a DCB-only approach removes this restriction, as no metal implant is left in the artery. Natural vasomotion supports positive remodelling and therefore late lumen enlargement can occur overtime<sup>8</sup>.




For a clinical paper relating to late lumen enlargement, please refer to reference 1.


## WHAT IS LATE LUMEN ENLARGEMENT?



Pre-treatment



Post-treatment



4-month follow-up

**Patient:** Female, 67 years

**Pre-treatment:** De novo stenosis of obtuse marginal branch

**Procedure:** **LESION PREPARATION**

Pre-dilation PTCA balloon (2.5 x 15 mm)

**LESION THERAPY**


DCB-Only SeQuent® Please (2.5 x 20 mm)

**Post-treatment:** Residual stenosis of <30%, slightly hazy in dilated stenosis but good TIMI flow

**Follow-up:** Perfect result without any restenosis demonstrating late lumen gain over time

Full case studies available on request.


# DCB-ONLY APPROACH CONSENSUS GUIDELINES<sup>9</sup>




# DCB-ONLY APPROACH CONSENSUS GUIDELINES

In comparison to stent implantation, the DCB-only strategy may take longer as adequate lesion preparation must be done gently and with patience. However, the patient can go home without any metalwork in their arteries and therefore has no risk of stent related complications which may require a second intervention.


## KEY POINTS:



Adequate lesion preparation is mandatory when using the DCB-only approach



The DCB is solely used as a drug delivery system



The DCB can be inflated again and used as a standard PTCA balloon if required



# DCB-ONLY APPROACH

## THE ANGIOPLASTY PROCESS

### STEP ONE

Mandatory lesion preparation using a standard PTCA balloon. The standard PTCA balloon is advanced to the lesion site.




### STEP TWO

The standard PTCA balloon is inflated, causing the lesion to crack.




### STEP THREE

The standard PTCA balloon is deflated and removed.




### STEP FOUR

The DCB is advanced to the lesion site as quickly as possible. As soon as the balloon enters the bloodstream, the drug, paclitaxel will become active.




### STEP FIVE

The DCB is inflated for 30 seconds at nominal pressure, paclitaxel with the help of its excipient iopromide, is instantly absorbed into the vessel wall.




### STEP SIX

The DCB is deflated and removed from the vessel, leaving nothing behind.





### STEP SEVEN

1 month DAPT regime begins<sup>8,9,10</sup>



### FOLLOW UP

If follow up is performed, positive remodelling and late lumen enlargement of the target vessel can be seen<sup>8</sup>



## SEQUENT® PLEASE NEO



### CLINICALLY PROVEN SECOND GENERATION DCB CATHETER WITH:

- Advanced crossing properties
- Improved pushability and deliverability
- Hydrophilic shaft coating
- Improved balloon profile and reduced balloon wall thickness
- Paclitaxel and Iopromide coating
- Effective drug release into the vessel wall
- Homogenous drug delivery

## SEQUENT® PLEASE NEO HANDLING TIPS

- Pre-dilate the target lesion with adequate lesion preparation
- The DCB is intended for a single dose transfer of paclitaxel
- The DCB can be used as a standard PTCA balloon following delivery of the drug, if required
- Sizing: 2 – 4 mm longer than pre-dilated area both distal and proximal to the lesion
- Balloon-vessel-diameter ratio: 1.0-1.0.
- Position the DCB at the lesion site as quickly as possible after entering the bloodstream
- Inflate once at nominal pressure for 30 seconds
- Do not bend, touch or wipe the balloon
- Do not expose the balloon to any fluids

# SEQUENT® PLEASE NEO ADVANTAGES

## NO UNNECESSARY STENT IMPLANTATION:

- No inflammation due to a foreign body implant
- No risk of stent related complications e.g. ISR or stent thrombosis
- No stent-related limitations for further treatment, if required
- No stent edge effect

## EFFICACY OF DCB:

- Enables positive remodelling<sup>1</sup>
- Keep natural vessel vasomotion
- 1 month DAPT<sup>1,2,3</sup>

# WHAT IS THE KEY MESSAGE?



Drug Coated Balloons will not completely replace drug eluting stents – drug eluting stents will continue to play a crucial role within coronary angioplasty.



DCBs can be used successfully in many indications outside of ISR and therefore are an alternative treatment option without stent-related disadvantages.



Patient-centered outcomes remain at the heart of clinical care and treatment options. Across Europe and the UK, the use of the DCB-only approach is increasing. Therefore, if a patient was eligible for the DCB-only approach and could avoid a metal implant, is it always necessary to implant one?

# SEQUENT® PLEASE NEO

## ORDERING INFORMATION

Balloon Diameter	Balloon Length						Nominal Pressure	Rated burst Pressure
	15 mm	20 mm	25 mm	30 mm	35 mm	40 mm		
2.0 mm	5023210	5023220	5023230	5023240	5023250	5023260	6 atm	14 atm
2.5 mm	5023212	5023222	5023232	5023242	5023252	5023262	6 atm	14 atm
3.0 mm	5023214	5023224	5023234	5023244	5023254	5023264	6 atm	14 atm
3.5 mm	5023216	5023226	5023236	5023246	5023256	5023266	6 atm	14 atm
4.0 mm	5023217	5023227	5023237	5023247	5023257	5023267	6 atm	14 atm

# SEQUENT® PLEASE NEO

## TECHNICAL DATA

Proximal shaft	1.9 F
Distal shaft	2.5 F
Usable length	145 cm
Balloon crossing profile	0.033" - 0.037"
Lesion entry profile	0.016"
Guiding catheter compatibility	5 F standard guiding catheter
Guidewire compatibility	0.014"
Rated burst pressure [RBP]	14 atm

# REFERENCES

1. Nordmann AJ, Briel M, Bucher HC. Mortality in randomized controlled trials comparing drug-eluting vs. bare metal stents in coronary artery disease: A meta-analysis. *Eur Heart J* 2006;27:2784–2814.
2. Camenzind E. Safety of drug-eluting stents: Insights from meta analysis. Presented at: European Society of Cardiology Annual Congress; September 2–6, 2006; Barcelona, Spain. (Hotline I; 707009)
3. Palmerini T, Benedetto U, Biondi-Zoccai G, Della Riva D, Bacchi-Reggiani L, Smits P, Vlachojannis G, Jensen L, Christiansen E, Berencsi K, Valgimigli M, Orlandi C, Petrou M, Rapezzi C and Stone G, 2015. Long-Term Safety of Drug-Eluting and Bare-Metal Stents. *Journal of the American College of Cardiology*, 65(23), pp.2496–2507.
4. Jeger R, Farah A, Ohlow M, Mangner N, Möbius-Winkler S, Leibundgut G, Weilenmann D, Wöhrle J, Richter S, Schreiber M, Mahfoud F, Linke A, Stephan F, Mueller C, Rickenbacher P, Coslovsky M, Gilgen N, Osswald S, Kaiser C, Scheller B, Buser P, Kühne M, Zellweger M, Sticherling C, Wein B, Twerenbold R, Fahrni G, Plicht B, Struck B, Önal I, Cremers B, Clever Y, Ewen S, Schirmer S, Böhm M, Wagner A, Lauer B, Stachel G, Höllriegel R, Winzer E, Rickli H, Ammann P, Haager P, Trachsel L, Joerg L, Nüssli D, Roelli H, Maeder M, Rohner F, Markovic S, Paliskyte R, Buckert D, Awad B, Erne P, Jamshidi P, Cuculi F, Kapos I, Toggweiler S, Riede F, Pörner T, Lenk K, Noutsias M, Surber R, Dannberg G, Franz M, Otto S, Zweiker R, Nlederl E, Perl S, Pieske B, Schmidt A, Luha O, Von Lewinski D, Krackhardt F, Kherad B, Jerichow T, Butter C, Neuss M, Tambor G, Hölschermann F, Bruch L, Winkler S, Lenz C, Seidel M, Keweloh B, Röttgen A, Bohl S, Wolf A and Hoffmann A, 2018. Drug-coated balloons for small coronary artery disease (BASKET-SMALL 2): an open-label randomised non-inferiority trial. *The Lancet*, 392(10150), pp.849–856.
5. Pleva L, Kukla P and Hlinomaz O, 2018. Treatment of coronary in-stent restenosis: a systematic review. *Journal of Geriatric Cardiology*, [online] 15(5), pp.173–184. Available at: <<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5895957/#>>.
6. Nice.org.uk. 2020. Medical Technologies Consultation | Sequent Please Balloon Catheter For In-Stent Coronary Restenosis | Guidance | NICE. [online] Available at: <<https://www.nice.org.uk/guidance/mtg1>>.
7. Escardio.org. 2020. ESC/EACTS Guidelines On Myocardial Revascularization. [online] Available at: <<https://www.escardio.org/Guidelines/Clinical-Practice-Guidelines/ESC-EACTS-Guidelines-in-Myocardial-Revascularisation-Guidelines-for>>.
8. Kleber F, Schulz A, Waliszewski M, Hauschild T, Böhm M, Dietz U, Cremers B, Scheller B and Clever Y, 2014. Local paclitaxel induces late lumen enlargement in coronary arteries after balloon angioplasty. *Clinical Research in Cardiology*, 104(3), pp.217–225.
9. Raban V, Jeger R, Simon Eccleshall, Wan Azman Wan Ahmad, Junbo Ge, Tudor C, Poerner, Eun-Seok Shin, Fernando Alfonso, Azeem Latib, Paul J. Ong, Tuomas T. Rissanen, Jorge Saucedo, Bruno Scheller, Franz X. Kleber, Drug-Coated Balloons for Coronary Artery Disease: Third Report of the International DCB Consensus Group, *JACC: Cardiovascular Interventions*, Volume 13, Issue 12, 2020, Pages 1391–1402
10. Corballis N, Wickramarachchi U, Vassiliou V and Eccleshall S, 2019. Duration of dual antiplatelet therapy in elective drug-coated balloon angioplasty. *Catheterization and Cardiovascular Interventions*.

**B. Braun Australia Pty Ltd**

Level 5, 7–9 Irvine Place, Bella Vista, NSW 2153 Australia

Tel.1800 251 705 | [info.au@bbraun.com](mailto:info.au@bbraun.com) | [www.bbraun.com.au](http://www.bbraun.com.au)