

TDS – Prontosan® Wound Gel | Prontosan® Wound Gel X

Cleansing and moistening of skin wounds and burns. For the prevention of biofilm.



- Removes & prevents biofilm
- Prevents infections
- Reduces healing time
- Painless dressing changes
- Well-known substances with low allergenic potential
- Compatible with commonly used wound dressings
- Can be used up to 8 weeks after first opening

While necrotic tissue and slough can be optimally removed by Prontosan® Wound Irrigation Solution, it is necessary to first remove heavy encrustations by keeping the wound moist for a longer period of time. For this task, Prontosan® Wound Gel and Gel X provide the optimal solution.

WHEN TO USE WHICH GEL

Prontosan® Wound Gel	Prontosan® Wound Gel X
30 ml pod	50 g tube
Fluid	Highly viscous
For smaller, deep and tunneling wounds, as well as for difficult to access wounds. For superficial and superficial partial thickness burns.	For large surface area wounds, e.g. <ul style="list-style-type: none"> • Ulcus cruris • Decubitus For deep partial thickness and full thickness burns.
If little quantity is required	When large quantities are required where Prontosan® Wound Gel is too fluid

STORAGE CONDITIONS

Store +5 to +30 °C

SHELF-LIFE

36 months

STABILITY



Can be used up to 8 weeks after first opening

TDS – Prontosan® Wound Gel | Prontosan® Wound Gel X

Cleansing and moistening of skin wounds and burns. For the prevention of biofilm.

HINTS AND TIPS

All wounds should in principle first be rinsed and cleansed with Prontosan® Wound Irrigation Solution, Prontosan® Gel and Gel X remain on the wound until the next dressing change. It therefore has a long acting effect.



INDICATIONS

Cleansing, moistening and decontamination of skin wounds and burns

- Traumatic wounds
- Postoperative wounds
- Chronic skin ulcers (e.g. venous, diabetic or pressure injuries)
- Thermal wounds
- Chemical wounds (acid- and alkali-induced)
- Radiation-induced wounds

Moistening of bandages and wound dressings



APPLICATION

The choice between Prontosan® Wound Gel and Gel X allows an optimal application in large surface area as well as deep wounds.

- 1 For large surface area wounds apply a 3-4 mm thick layer of the Prontosan® Wound Gel X and cover with a secondary dressing.
- 2 For the application in deep or tunneling wounds, wound cavities and difficult to access areas, apply a 3-5 mm layer of Prontosan® Wound Gel and cover with a secondary dressing.

Prontosan® Wound Gel



	Size	Pcs / Pack	Article number
Prontosan® Wound Gel	30 ml	20	400510
Prontosan® Wound Gel X	50 g	20	400523

Composition: 0.1 % Polyaminopropyl Biguanide (Polihexanide), Betaine Surfactant, Purified Water, Glycerol, Hydroxyethylcellulose