

# ECOTAINER®

## IRRIGATION CONTAINER

### Solutions well in hand

Whenever you need larger volumes of irrigation solution, you can rely on Ecotainer®.

The container is made of polypropylene - a transparent and robust material suitable to meet the highest demands.

The moulded shape of Ecotainer® offers a comfortable grip for all hands regardless of size, even under wet conditions.



The colour coded screw cap is easy to open. The white sealing ring gives visible assurity that the container is sealed and tamper-proof.

- Sodium Chloride 0.9% for irrigation
- Water for irrigation



Ecotainer® provides excellent non-drip pouring characteristics which allow irrigation solutions to be applied cleanly and accurately.

A graduation scale supports volume measurement.

Fluid	Size	Box Qty	B. Braun Code	CH2 Code	Symbion
<span style="display: inline-block; width: 15px; height: 15px; border: 1px solid #0070C0; margin-right: 5px;"></span> Water for irrigation	250ml	12	3637007	2470831	665320
	500ml	10	3521380	2069618	130206
	1000ml	6	3521390	2069600	130184
<span style="display: inline-block; width: 15px; height: 15px; background-color: #0070C0; margin-right: 5px;"></span> Sodium Chloride 0.9% for irrigation	250ml	12	3637006	2470823	665339
	500ml	10	3521360	2069642	130176
	1000ml	6	3521370	2069634	130168