

Aesculap Endoscopic Technology Suction and Irrigation System



Suction/Irrigation Device

Saug- Spülsystem



Features

- Colour coded buttons
- Transparent valve
- Ergonomically designed handle
- Double tubing with spike connection for irrigation bags
- Blunt tip probe
- Distal holes – maximum suction performance
- Quick and easy assembling/ disassembling of the probe

Eigenschaften

- Farbcodierte Saug/Spülknöpfe
- Transparentes Ventil
- Ergonomisches Handstück
- Doppeltes Schlauchsystem mit Dornen zum Anschluss an Spülbeutel
- Saug/Spülrohr mit stumpfer Spitze
- Distale Bohrungen – maximale Saugleistung
- Schnelle und einfache Montage/ Demontage des Saug/Spülrohrs

Cat. No. Art. Nr.	Description Beschreibung	Diameter Durchmesser	Length Länge	Quantity Menge
Complete package including probe / Set komplett mit Saug-Spül-Rohr				
PG036SU	Disposable Suction-Irrigation probe with valve and double spike connection for irrigation bags, sterile Einmal Saug-Spül-Rohr und Handstück mit doppeltem Dorn für Spülbeutel, steril	5 mm	330 mm	10 pieces 10 Stück
Valves only – for use with reusable probes / Handstück alleine – zur Verwendung mit wiederverwendbarem Saug-Spül-Rohr				
PG037SU	Disposable valve with double spike connection for irrigation bags, sterile Einmal Handstück mit doppeltem Dorn für Spülbeutel, steril			10 pieces 10 Stück
Reusable Probes / Saug-Spül-Rohre, wiederverwendbar				
PG038R	Suction-Irrigation probe with distal holes Saug-Spül-Rohr mit distalen Bohrungen	5 mm	330 mm	1 piece 1 Stück
PG039R	Suction-Irrigation probe with distal holes Saug-Spül-Rohr mit distalen Bohrungen	10 mm	330 mm	1 piece 1 Stück



AESCULAP®

B | BRAUN
SHARING EXPERTISE

Aesculap AG & Co. KG

Am Aesculap-Platz
78532 Tuttlingen
Germany

Phone +49 7461 95-0
Fax +49 7461 95-2600

www.aesculap.de

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogues and pricelists and to take legal actions.