

Introcan Safety[®] 3 in radiology

Improving safety and efficiency in your daily work.



User and Patient Benefits

- The integrated Stabilization Platform is designed to minimize catheter movement and helps reducing irritation and restarts.¹⁹
- B. Braun's Introcan Safety[®] 3 Closed IV Catheters 24–18 gauge are suitable for use with power injectors set to a maximum pressure of 300 psi.¹⁵
- 22G accommodates to a commonly used flow rate from 3 mL/s to 5 mL/s.¹⁶ The smallest gauge and the shortest length of cannula that will accommodate the prescribed therapy should be selected according to nursing recommendations. A small gauge cannula results in less trauma to the vessel and allows adequate blood flow around the catheter wall.¹⁸

