

STERICAN® NEEDLES

SINGLE USE HYPODERMIC NEEDLES

Size	Needle Length	Outer Ø	Box Qty	Code
Sterican® Conventional Hypodermic Needles				
20G	1 1/2"	40 mm	0.90 mm	100 4657519
21G	1 1/2"	40 mm	0.80 mm	100 4657527
22G	1 1/4"	30 mm	0.70 mm	100 4657624
23G	1 1/4"	30 mm	0.60 mm	100 4657640
23G	1"	25 mm	0.60 mm	100 4657667
24G	1"	25 mm	0.55 mm	100 4657675
26G	1"	25 mm	0.45 mm	100 4657683
27G	3/4"	20 mm	0.40 mm	100 4657705
Sterican® Insulin Needles				
26G	1/2"	12 mm	0.45 mm	100 4665457
27G	1/2"	12 mm	0.40 mm	100 4665406
30G	1/2"	12 mm	0.30 mm	100 4656300
Sterican® Heparin, Tuberculin				
25G	5/8"	16 mm	0.50 mm	100 4657853
Sterican® for Blood Sampling				
19G	1 1/2"	40 mm	1.10 mm	100 4657799
19G	1 1/4"	30 mm	1.10 mm	100 4665317
20G	1"	25 mm	0.90 mm	100 4657500
21G	1"	25 mm	0.80 mm	100 4657543
22G	1 1/2"	40 mm	0.70 mm	100 4660021
Sterican® for Phlebotomy/Blood Sampling				
14G	3 1/8"	80 mm	2.10 mm	100 4665473
18G	1 1/2"	40 mm	1.20 mm	100 4665120
18G	2"	50 mm	1.20 mm	100 4667123
19G	1 1/2"	40 mm	1.10 mm	100 4665112
20G	2"	50 mm	0.90 mm	100 4667093
Sterican® Deep Intramuscular, Long Bevel				
20G	2 3/4"	70 mm	0.90 mm	100 4665791
21G	3 1/8"	80 mm	0.80 mm	100 4665465
Sterican® Intramuscular, Long Bevel				
19G	2"	50 mm	0.10 mm	100 4665511
21G	2"	50 mm	0.80 mm	100 4665503
22G	2"	50 mm	0.70 mm	100 4650018
Sterican® Dental				
25G	1"	25 mm	0.50 mm	100 9186158
25G	1 1/2"	40 mm	0.50 mm	100 9186166
27G	1"	25 mm	0.40 mm	100 9180117
27G	1"	25 mm	0.40 mm	100 9186174
Sterican® for Neural Therapy				
21G	4 3/4"	120 mm	0.80 mm	100 4665643
23G	2 3/8"	60 mm	0.60 mm	100 4665600



FEATURES:

- Thin-walled needles for minimal pain upon skin puncture
- Silicone coating makes needle surface extremely smooth to minimise discomfort
- Transparent Luer-Lock plastic hub for easy visibility
- Conforms to ISO 7864, ISO 6009 and DIN 13097

STERICAN® MIX BLUNT FILL NEEDLES FOR DRUG ADMIXTURE ONLY



- *40° and **90° blunt fill needles offer clinicians flexibility
- Blunt and semi blunt needles minimises the risk of needlestick injuries
- Thin walled cannula allows for faster medication draw up
- *Pink coloured hub and protective cap for blunt fill needles
- **Red coloured hub and protective cap for semi blunt needles

Size	Needle Length	Outer Ø	Carton	Code
* 18G	1 1/2"	40 mm	1000	4038088-03
** 18G	1 1/2"	40 mm	1000	4550400-03

SYRINGES

INJEKT® TWO-PIECE STANDARD SYRINGES



- Green plunger for easy recognition of filling volume
- Black graduations for ideal contrast and readability
- High barrel transparency simplifies inspection of the content
- Polypropylene syringe barrel and plunger
- Latex and DEHP Free

Description	Size	Box Qty	Code
Luer Slip	1 mL	100	9166017V
Luer Slip	2 mL	100	4606027V
Luer Slip	5 mL	100	4606051V
Luer Slip	10 mL	100	4606108V
Luer Slip	20 mL	100	4606205V
Luer Lock	2 mL	100	4606701V
Luer Lock	5 mL	100	4606710V
Luer Lock	10 mL	100	4606728V
Luer Lock	20 mL	100	4606736V

OMNIFIX® THREE-PIECE STANDARD SYRINGES



- Plunger stopper with double sealing ring
- Easily noticeable backstop to prevent unintentional withdrawal of plunger from barrel
- Colour-coded dispenser boxes for quick identification by size
- Low stiction properties
- Polypropylene barrel and plunger
- Latex and DEHP free

Description	Size	Box Qty	Code
Luer Slip	1 mL	100	9161406V
Luer Slip	3 mL	100	4616025V
Luer Slip	5 mL	100	4616057V
Luer Slip	10 mL	100	4616103V
Luer Slip	20 mL	100	4616200V
Luer Slip	30 mL	100	4616308F
Luer Slip	50 mL	100	4616502F
Luer Slip	50/60 mL	100	4613503F
Luer Lock	3 mL	100	4617022V
Luer Lock	5 mL	100	4617053V
Luer Lock	10 mL	100	4617100V
Luer Lock	20 mL	100	4617207V
Luer Lock	30 mL	100	4617304F
Luer Lock	50/60 mL	100	4617509F

OMNICAN® N INSULIN SYRINGES



- Highly transparent barrel with high contrast black graduation for easy reading and exact dosing
- Highly transparent barrel to make sure that air bubbles are easily identified
- Integrated thin-wall needle which allows greater comfort when injecting

Description	Size	Box Qty	Code
Omnican® N 50	30G x 8 mm	100	FT9172210S
Omnican® N 100	30G x 8 mm		FT9172220S
Omnican® N 100	30G x 12 mm		FT9172230S
Omnican® N 40	30G x 12 mm		FT9172240S